



The Code ORANGE XL sharps container offers the following features:-

MANUFACTURED, TESTED AND CERTIFIED TO BS7320:1990

UN TESTED AND CERTIFIED

ORANGE COLOUR CODED CLOSURE TO MEET THE RECOMMENDATIONS ON WASTE SEGREGATION AS DETAILED IN THE DEPARTMENT OF HEALTH'S TECHNICAL MEMORANDUM HTM-07-01 – SAFE MANAGEMENT OF HEALTH CARE WASTE.

EASY TO OPERATE SLIDING LID

SECURE CLOSURE LOCK

DISTINCTIVE COLOUR

TOUGH AND DURABLE

SUPPLIED ASSEMBLED

ONE PIECE BASE INCORPORATING HANDLES

COST EFFECTIVE

PURPOSE DESIGNED KICK TROLLEY

Sharpak XL sharps container

Following the publication of the Department of Health Technical Memorandum HTM-07-01 titled 'The Safe Management of Healthcare Waste' Amcor Flexibles are able to offer a range of colour coded sharps containers to assist in the segregation of waste as recommended in the new guide.

The SHARPAK XL is a specialist container aimed at theatre sharps disposal, in particular, where minimal invasive surgery procedures involving the use of long disposable instruments. Its large aperture combined with its overall extended body length facilitates the disposal of the very longest and bulkiest endoscopy instruments, trocars, chest drains, pipettes and other sharps.

Code ORANGE Sharpak sharps containers can be disposed of by incineration or other approved method.

Please note:- For disposal of sharps contaminated with cytotoxic or cytostatic medicines please refer to the Amcor Flexibles SPS data sheet for the Code PURPLE range of Sharpak HYDRI sharps containers.



The SHARPAK XL is approved to the British Standard BS7320:1990 and tested and certified in accordance with the provisions on the Transportation of Dangerous Goods, Chapter 9, UN3291 Clinical Waste Unspecified, NOS.

SHARPAK™
HEALTHCARE




SHARPAK® PRODUCT RANGE

PRODUCT	MAX CAPACITY	CAPACITY TO FILL LINE	Re-order code	HEIGHT	WIDTH	LENGTH
Sharpak XL	26.3 litres	19.2 litres	SU06520	230 mm	190 mm	660 mm



Dimensions in mm

Actual aperture sizes when printed at 100%



For further information on any of these Sharpak sharps containers see the product specific data sheets.

For details regarding SHARPAK distributors and world wide availability, please contact:
Amcor Flexibles/SPS Laboratoires Winterbourne Road Stoke Gifford BRISTOL BS34 6PT
Tel +44 (0)117 983 6000 Fax +44 (0)117 983 6001
www.sharpakforsharps.co.uk

