



Sharpak COMPACT, COMPACT 'O', 6, 18 PLUS and 18 Eco sharps containers

Following the publication of the Department of Health Technical Memorandum HTM-07-01 titled 'The Safe Management of Healthcare Waste' Amcor Flexibles are able to offer a range of colour coded sharps containers to assist in the segregation of waste as recommended in the new guide.

All Sharpak Code YELLOW sharps containers feature yellow coloured closures and are recommended for the segregation, collection and disposal of partially discharged sharps waste which are NOT contaminated with cytotoxic or cytostatic medicines.

Code YELLOW Sharpak sharps containers MUST only be disposed of by incineration.

Please note:- For disposal of sharps contaminated with cytotoxic or cytostatic medicines please refer to the Amcor Flexibles SPS data sheet for the Code PURPLE range of Sharpak HYDRI sharps containers.

All Sharpak Code YELLOW sharps containers are approved to BS7320:1990 and tested and certified in accordance with the provisions of the United Nations recommendations on the Transportation of Dangerous Goods, Chapter 9, UN3291 Clinical Waste, Unspecified, NOS.

The Code YELLOW range of sharps containers offers the following features:-

MANUFACTURED, TESTED AND CERTIFIED TO BS7320:1990

UN TYPE TESTED

YELLOW COLOUR CODED CLOSURES TO MEET THE RECOMMENDATIONS ON WASTE SEGREGATION AS DETAILED IN THE DEPARTMENT OF HEALTH'S TECHNICAL MEMORANDUM HTM-07-01 – SAFE MANAGEMENT OF HEALTHCARE WASTE

SECURE CLOSURE WITH LOCKING FACILITY

VISILOCK INDICATOR IS FITTED ON THE SHARPAK COMPACT, COMPACT 'O' AND SHARPAK 18 PLUS

DISTINCTIVE COLOURING AND CLEAR LABELLING

SEQUENTIAL NUMBERING FOR TRACEABILITY

TOUGH AND DURABLE

SUPPLIED ASSEMBLED

COST EFFECTIVE

PURPOSE DESIGNED WALL/TROLLEY BRACKETS AND COMPACT STAND

SUPPLIED WITH SHARPAK HYDRI MATS TO ABSORB ANY SMALL AMOUNTS OF FREE LIQUID IN THE EVENT OF AN ACCIDENTAL SPILLAGE IN THE CONTAINER

INSULIN PEN NEEDLE REMOVING FEATURE ON SHARPAK 18 ECO




SHARPAKTM
HEALTHCARE



VisiLockTM
PATENTED

SHARPAK® PRODUCT RANGE

PRODUCT	MAX CAPACITY	CAPACITY TO FILL LINE	Re-order code	HEIGHT	WIDTH	DEPTH
Sharpak COMPACT	0.5 litres	0.39 litres	SU06230	186 mm	67 mm	67 mm
Sharpak COMPACT 'O'	0.5 litres	0.39 litres	SU06232	186 mm	67 mm	67 mm
Sharpak 6	0.6 litres	0.48 litres	SU07430	200 mm	102 mm	46 mm
Sharpak 18 PLUS	1.58 litres	1.27 litres	SU06430	210 mm	104 mm	104 mm
Sharpak 18 Eco	1.53 litres	1.22 litres	SU06432	190 mm	104 mm	104 mm



Actual aperture sizes when printed at 100%

Dimensions in mm

For further information on any of these Sharpak sharps containers see the product specific data sheets.

For details regarding SHARPAK distributors and world wide availability, please contact:
Amcor Flexibles/SPS Laboratoires Winterbourne Road Stoke Gifford BRISTOL BS34 6PT
Tel +44 (0)117 983 6000 Fax +44 (0)117 983 6001
www.sharpakforsharps.co.uk

