

Sharpak 18 Eco sharps container

As an alternative to the standard Sharpak 18 sharps container, the 1.53 litre SHARPAK 18® Eco has been developed to offer a sharps container that is not only compact but extremely versatile. Made from a robust plastic, the SHARPAK 18 Eco is suitable for the vast majority of sharps disposal needs without occupying valuable space. Pre-assembled for convenience, this container offers all the qualities normally associated with a UN/BS SHARPAK product, including an increased aperture size which is capable of receiving larger items such as transfer devices. Also included are needle, vacutainer and insulin pen needle removing features, and a strong integral handle to suit wall and trolley mounting.

The SHARPAK 18 Eco offers the following features:

MANUFACTURED, TESTED AND CERTIFIED TO BS7320:1990
AND FRENCH STANDARD NF X 30-500

UN TESTED AND CERTIFIED

SUPPLIED ASSEMBLED

EASY TO OPERATE CLOSURE AND FULLY LOCKABLE LID

DISTINCTIVE COLOURING AND CLEAR LABELLING

TOUGH AND DURABLE

SEQUENTIAL NUMBERING FOR TRACEABILITY

COST EFFECTIVE

PURPOSE DESIGNED WALL AND TROLLEY BRACKET


INCREASED APERTURE CAPABLE OF ACCEPTING
LARGER ITEMS INCLUDING TRANSFER DEVICES

INSULIN PEN NEEDLE REMOVING FEATURE


SUITABLE FOR USE WITH THE SHARPAK PROCEDURE TRAY

The SHARPAK 18 Eco is approved to the British Standard BS7320:1990 and tested and certified in accordance with the provisions of the United Nations recommendations on the Transportation of Dangerous Goods, Chapter 9, UN3291 Clinical Waste, Unspecified, NOS.






TOP VIEW




Actual aperture size
Dimensions in mm



Insulin pen needle
removing feature



BACK VIEW



SIDE VIEW

SHARPAK® 18 ECO INFORMATION

Height	190 mm
Width	104 mm
Aperture	49.5 mm dia.
Max Capacity	1.53 litres
Capacity to fill line	1.22 litres
Sharpak 18 Eco case re-order code	SU06496
Quantity per case	24
Wall and trolley bracket	SM06500



OPEN



PUSH HALF WAY TO TEMPORARY CLOSE



PUSH DOWN FULLY TO LOCK



REMOVAL OF
DISPOSABLE
INSULIN PEN
NEEDLE

For details regarding SHARPAK distributors and world wide availability, please contact: Amcor Flexibles/SPS Laboratoires
Winterbourne Road Stoke Gifford BRISTOL BS34 6PT
Tel +44 (0)117 983 6003
Fax +44 (0)117 983 6001
E mail david.scarrow@amcor-flexibles.com
www.sharpakforsharps.co.uk