

Universal mounting bracket



WALL MOUNTED



STAND/TROLLEY MOUNT


The universal accessory pack contains all the necessary parts to allow the SHARPAK® 60, 120 and 220 PLUS containers to be safely mounted either on a wall or fixed to a trolley or equipment stand.

BRACKET INFORMATION

Pack Quantity	1
Re-order code	SM05000



Sharpak kick trolley



The Kick Trolley from Sharpak Healthcare has been designed to accept the SHARPAK® 120 and 220 PLUS sharps containers allowing them to be safely located at a sensible work height with the added benefit of being able to be moved around the department when and where needed. One wheel is fitted with a lock.

KICK TROLLEY INFORMATION

Pack Quantity	1
Re-order code	SM16046



Sharpak XL trolley

Allows use of the Sharpak XL container in its vertical position and can be easily moved around within the procedure area.

NEW XL TROLLEY INFORMATION

Pack Quantity	1
Re-order code	SM16047

For details regarding SHARPAK distributors and world wide availability, please contact: Amcor Flexibles/SPS Laboratoires
Winterbourne Road Stoke Gifford BRISTOL BS34 6PT
Tel +44 (0)117 983 6003
Fax +44 (0)117 983 6001
E mail david.scarrow@amcor-flexibles.com
www.sharpakforsharps.co.uk