



## MAGPIE™

The MAGPIE waste container provides a simple, safe and cost effective way to dispose of 'empty' aerosols and broken glass waste in hospitals, clinics, laboratories and kitchens.

The sturdy corrugated board container is supplied with a polyethylene coated solid board insert to prevent perforation and leakage. The container is clearly marked with assembly and closure instructions to prevent misuse.

The MAGPIE container and insert are supplied flat packed to save storage space and reduce transport cost.

The MAGPIE offers the following features:-

---

SUPPLIED FLAT PACKED

---

STRONG CONSTRUCTION

---

SIMPLE TO ASSEMBLE AND CLOSE

---

SUPPLIED WITH A POLYETHYLENE COATED SOLID BOARD INSERT FOR ADDED STRENGTH

---

INCORPORATED CARRY HANDLE

---

STACKABLE

---

FACILITY TO ALLOW THE LID TO BE SECURED USING A STANDARD PLASTIC SECURITY TIE

---

### MAGPIE INFORMATION

---

Max Capacity 24 litres

---

Height 450 mm

---

Width 275 mm

---

Depth 275 mm

---

SHARPAK 2 re-order code CW04165

---

Pack quantity 10

For details regarding SHARPAK distributors and world wide availability, please contact: Amcor Flexibles/SPS Laboratoires Winterbourne Road Stoke Gifford BRISTOL BS34 6PT  
Tel +44 (0)117 983 6003  
Fax +44 (0)117 983 6001  
E mail david.scarrow@amcor-flexibles.com  
www.sharpakforsharps.co.uk

