

In response to specific requests from our customers the Sharpak PHARMApak containers feature green colour coded closures as well as green labels to aid container identification during waste segregation and disposal.

CODE GREEN



## Sharpak 60, 120 and 220 PHARMApak containers with VisiLOCK.

The Sharpak 60, 120 and 220 PHARMApak containers, as well as being approved for the collection and disposal of used sharps, have now also been UN type tested and approved for the collection of waste medicines in their original packagings.

The UN type testing was carried out using test fills consisting of a mixture of medicines including tablets, capsules, creams, gels, powders, pastes, liquid unit dose drugs, ampoules/vials, bottles of liquid and syrups in their original packagings and needles and syringes. It is the responsibility of the consignor to record the correct UN transport classification number or European Waste Code number on to the container label e.g. UN 3249, UN 1851 and/or EWC code 18 01 09. To obtain the correct UN classification or European Waste Code (EWC) for the waste in question, consult your Dangerous Goods Safety Advisor (DGSA).

These Code GREEN Sharpak PHARMApak containers MUST NOT be used for the collection and disposal of any cytotoxic or cytostatic waste. For dealing with this type of healthcare waste please refer to the product data sheet which covers the Code PURPLE Sharpak HYDRI range of sharps containers.

The Sharpak 60, 120 and 220 PHARMApaks offer the following features:-

MANUFACTURED AND PERFORMANCE TESTED TO BS7320:1990

UN TYPE TESTED FOR SHARPS WASTE

UN TYPE TESTED FOR WASTE MEDICINES

GREEN COLOUR CODED CLOSURES AND LABELS TO AID CONTAINER IDENTIFICATION DURING WASTE SEGREGATION AND DISPOSAL

EASY TO OPERATE ROTATING OPEN/CLOSE AND FINAL LOCK CLOSURE

SECURE CLOSURE WITH VISILOCK INDICATOR

DISTINCTIVE COLOURING AND CLEAR LABELLING

SEQUENTIAL NUMBERING FOR TRACEABILITY

TOUGH AND DURABLE

QUICK AND EASY TO ASSEMBLE

COST EFFECTIVE

PURPOSE DESIGNED WALL/TROLLEY BRACKETS AND KICK TROLLEY

SUPPLIED WITH SHARPAK HYDRI MATS TO ABSORB ANY SMALL AMOUNTS OF FREE LIQUID IN THE EVENT OF AN ACCIDENTAL SPILLAGE IN THE CONTAINER

**SHARPAK**<sup>TM</sup>  
HEALTHCARE



**VisiLOCK**<sup>TM</sup>  
PATENTED







In response to specific requests from our customers the Sharpak PHARMApak containers feature green colour coded closures as well as green labels to aid container identification during waste segregation and disposal.

**SHARPAK® PHARMApak INFORMATION**

PRODUCT	MAX CAPACITY	TESTED CAPACITY	Re-order code	Quantity Per Case	HEIGHT	WIDTH	DEPTH
PHARMApak 60	6 litres	2kg	PH02000	25	299 mm	220 mm	234 mm
PHARMApak 120	12 litres	3kg	PH02200	15	325 mm	239 mm	268 mm
PHARMApak 220	22 litres	5kg	PH02400	8	383 mm	297 mm	327 mm



*Dimensions in mm*

Actual aperture sizes when printed at 100%

For details regarding SHARPAK distributors and world wide availability, please contact: Amcor Flexibles/SPS Laboratoires  
Winterbourne Road Stoke Gifford BRISTOL BS34 6PT  
Tel +44 (0)117 983 6000 Fax +44 (0)117 983 6001  
[www.sharpakforsharps.co.uk](http://www.sharpakforsharps.co.uk)

