



## Sharpak Compact sharps containers with **VisiLOCK™** PATENTED

The SHARPAK® COMPACT sharps container has been specially developed to meet the growing needs of clinical staff involved in giving injections as well as for those who have to self administer medication by means of a disposable needle and syringe.

It is not always practical to have access to a safe means of disposing of used sharps and this container provides a purpose designed product which is safe, simple to use, portable and compact in size.

The SHARPAK COMPACT is approved to the British Standard BS7320:1990, French Standard NF X 30-500 and tested and certified in accordance with the provisions of the United Nations recommendations on the Transportation of Dangerous Goods, Chapter 9, UN3291 Clinical Waste, Unspecified, NOS.

The SHARPAK COMPACT offers the following features:

MANUFACTURED, TESTED AND CERTIFIED TO  
BS7320:1990 AND NF X 30-500

UN TESTED AND CERTIFIED

SUPPLIED ASSEMBLED

EASY TO OPERATE ROTATING OPEN/CLOSE/LOCK  
CLOSURE

APERTURE EQUIPPED WITH NON RETURN PETAL

SECURE CLOSURE LOCK

FEATURING VISILOCK INDICATOR SYSTEM

DISTINCTIVE COLOURING AND CLEAR LABELLING

TOUGH AND DURABLE

SEQUENTIAL NUMBERING FOR TRACEABILITY

COST EFFECTIVE

PURPOSE DESIGNED STAND AVAILABLE

ACCEPTS UP TO 5ML SYRINGES


SUITABLE FOR USE WITH THE SHARPAK PROCEDURE TRAY

ALSO AVAILABLE WITH AN OPEN APERTURE  
(NO NON-RETURN PETALS) COMPACT 'O'



There is also a demand for a sharps container with all the same safety and performance features as the SHARPAK COMPACT sharps container but where the user wishes to be more discreet. To meet this need the SHARPAK NX5 is available. The 'only' difference

between the two products is that the SHARPAK NX5 has been manufactured with a grey base and black closure and the labelling has been kept as small and as simple as possible.




Compact stand



OPEN



1ST POSITION 'HALF OPEN'



2ND POSITION 'CLOSED'

**SHARPAK® COMPACT INFORMATION**

|  |             |
|--|-------------|
| Max Capacity                           | 0.5 litres  |
| Capacity to fill line                  | 0.39 litres |
| SHARPAK COMPACT case re-order code     | SU06201     |
| SHARPAK COMPACT quantity per case      | 24          |
| SHARPAK COMPACT 'O' case re-order code | SU06202     |
| SHARPAK COMPACT 'O' quantity per case  | 24          |
| SHARPAK NX5 case order no              | SU06204     |
| SHARPAK NX5 quantity per case          | 24          |
| COMPACT/NX5 stand (pack quantity 1)    | SM07500     |

For details regarding SHARPAK distributors and world wide availability, please contact: Amcor Flexibles/SPS Laboratoires Winterbourne Road Stoke Gifford BRISTOL BS34 6PT  
Tel +44 (0)117 983 6000  
Fax +44 (0)117 983 6001  
E mail david.scarrow@amcor-flexibles.com  
www.sharpakforsharps.co.uk