



Sharpak procedure trays

SHARPAK® Procedure Trays have been specifically designed for use with the SHARPAK COMPACT, 18 & 35 'PLUS' sharps containers.

These practical devices provide clinical staff the means of safe handling and disposal of sharps when dealing with patients at the bedside, ensuring safe transportation of used sharps from the point of use, encouraging staff to adopt good practice.

The Sharpak range of Procedure Trays incorporate the following features:

STEAM STERILISABLE UP TO 134°C

DISTINCTIVE COLOUR ENSURES STAFF HAVE CONTINUED AWARENESS OF SHARPS HAZARD


LARGE TRAY AREA CAPABLE OF HOLDING THE NECESSARY MATERIALS REQUIRED TO CARRY OUT THE PROCEDURE

DEEP SIDES REDUCE THE RISK OF SPILLAGE

LOW UNIT COST


PACKED IN 5 UNITS PER CARTON






**SHARPAK 35 'PLUS' SHAPED
PROCEDURE TRAY**

| | |
|-----------------------------------|--------------|
| Depth | 60 mm |
| Size | 325 x 310 mm |
| Procedure Tray case re-order code | SU03501 |
| Quantity per carton | 5 |



**SHARPAK 18 & COMPACT SHAPED
PROCEDURE TRAY**

| | |
|-----------------------------------|--------------|
| Depth | 60 mm |
| Size | 310 x 285 mm |
| Procedure Tray case re-order code | SU06211 |
| Quantity per carton | 5 |



**SHARPAK 18 & COMPACT RECTANGULAR
PROCEDURE TRAY**

| | |
|-----------------------------------|--------------|
| Depth | 50 mm |
| Size | 208 x 308 mm |
| Procedure Tray case re-order code | SU06210 |
| Quantity per carton | 5 |



For details regarding SHARPAK distributors and world wide availability, please contact: Amcor Flexibles/SPS Laboratoires
Winterbourne Road Stoke Gifford BRISTOL BS34 6PT
Tel +44 (0)117 983 6003
Fax +44 (0)117 983 6001
E mail david.scarrow@amcor-flexibles.com
www.sharpakforsharps.co.uk